



Changzhou Meditech Technology Co., Ltd.

 Address: 9th floor, INNO LASER Building, Science & Education Town
Wujin District, Changzhou, China

 Tel: +86-519-85855955

 Fax: +86-519-85855955

 Website: www.czmeditech.com



MEDITECH
TECHNOLOGY
EXTERNAL FIXATION SERIES

2025

 www.czmeditech.com

COMPANY PROFILE

Changzhou Meditech Technology Co., Ltd., located in Changzhou Science & Education Town, Jiangsu province, is specialized in manufacturing orthopedic implants and instruments.

After more than 10 years research and development, We have 8 main product series and they are spinal system, intramedullary nail system, trauma plate and screw system, locking plate and screw system, CMF system, external fixator system, medical power tool system and general surgical instruments system.

With CE and ISO 13485 certificates, our products have been exported to USA, Germany, Argentina, Chile, Mexico, India, Thailand, Malaysia, Turkey, Egypt, South Africa, Ivory Coast, etc.

In the principle of "quality first, service first, R & D first, innovation first", the company wins an excellent reputation both in domestic and abroad. The company takes benefiting patients as its permanent goal and makes unremitting efforts for human health.



CONTENTS


1

Orthofix Type External Fixator System




11

Ilizarov Type External Fixator System



24

Hoffmann Type External Fixator System



38

Synthes AO Type External Fixation System




49

Common External Fixator System



59

Accessories of External Fixator System



Mini Fragment External Fixator-Type A

REF: 6100-0101

Matching bone screws: $\Phi 2 \times 70\text{mm}$ 4 pcs

Matching instruments: 3mm hex wrench, 2mm screwdriver

Length: 100mm




Mini Fragment External Fixator-Type B

REF: 6100-0102

Matching bone screws: $\Phi 2 \times 70\text{mm}$ 4 pcs

Matching instruments: 3mm hex wrench, 2mm screwdriver

Length: 90mm



Mini Fragment External Fixator-Type E

REF: 6100-0105

Matching bone screws: $\Phi 2 \times 70\text{mm}$ 4 pcs

Matching instruments: 3mm hex wrench, 2mm screwdriver

Length: 60mm



Mini Fragment External Fixator-Type F

REF: 6100-0106

Matching bone screws: $\Phi 2 \times 70\text{mm}$ 4 pcs

Matching instruments: 3mm hex wrench, 2mm screwdriver

Length: 70mm



Mini Fragment External Fixator-Type C

REF: 6100-0103

Matching bone screws: $\Phi 2.5 \times 50\text{mm}$ 4 pcs

Matching instruments: 3mm hex wrench, 2mm screwdriver

Length: 74mm



Mini Fragment External Fixator-Type D

REF: 6100-0104

Matching bone screws: $\Phi 2 \times 50\text{mm}$ 4 pcs

Matching instruments: 3mm hex wrench, 2mm screwdriver

Length: 45mm



Mini Fragment External Fixator-Type G

REF: 6100-0107

Matching bone screws: $\Phi 2 \times 70\text{mm}$ 4 pcs

Matching instruments: 3mm hex wrench, 2mm screwdriver

Length: 90mm



Mini Fragment External Fixator-Type H

REF: 6100-0108

Matching bone screws: $\Phi 2 \times 70\text{mm}$ 4 pcs

Matching instruments: 3mm hex wrench, 2mm screwdriver

Length: 70mm




Wrist External Fixator-Type A

REF: 6100-0201
Matching bone screws: $\Phi 3.5 \times 80\text{mm}$ 4 pcs
Matching instruments: 3.5mm screwdriver, 4mm hex wrench



Wrist External Fixator-Type B

REF: 6100-0202
Matching bone screws: $\Phi 2.5 \times 60\text{mm}$ 2 pcs, $\Phi 3.5 \times 80\text{mm}$ 2 pcs
Matching instruments: 10 mm open spanner, 2.5&3.5mm screwdriver, 4mm hex wrench




Wrist External Fixator-Type D

REF: 6100-0204
Matching bone screws: $\Phi 2.5 \times 60\text{mm}$ 2 pcs, $\Phi 3.5 \times 80\text{mm}$ 2 pcs
Matching instruments: 3.5mm screwdriver, 4mm hex wrench



Wrist External Fixator-Type E

REF: 6100-0205
Matching bone screws: $\Phi 2.5 \times 60\text{mm}$ 2 pcs, $\Phi 3.5 \times 80\text{mm}$ 2 pcs
Matching instruments: 3.5mm screwdriver, 4mm hex wrench



Wrist External Fixator-Type F

REF: 6100-0206
Matching bone screws: $\Phi 3.5 \times 80\text{mm}$ 4 pcs
Matching instruments: 3.5mm screwdriver, 5mm hex wrench, 10mm open spanner



Wrist External Fixator (Radiolucent Type)

REF: 6100-0207
Matching bone screws: $\Phi 2.5 \times 60\text{mm}$ 2 pcs, $\Phi 3.5 \times 80\text{mm}$ 2 pcs
Matching instruments: 3.5mm screwdriver, 4mm hex wrench



Wrist External Fixator (Radiolucent Type)

REF: 6100-0208
Matching bone screws: $\Phi 2.5 \times 60\text{mm}$ 2 pcs, $\Phi 3.5 \times 80\text{mm}$ 2 pcs
Matching instruments: 3.5mm screwdriver, 5mm hex wrench



Dynamic Axial External Fixator(Straight Type)

XL

REF: 6100-0301

Matching bone screws: $\Phi 6 \times 150\text{mm}$ 4 pcs

Matching instruments: 6mm hex wrench, 6mm screwdriver

L

REF: 6100-0302

Matching bone screws: $\Phi 6 \times 150\text{mm}$ 4 pcs

Matching instruments: 6mm hex wrench, 6mm screwdriver

M

REF: 6100-0303

Matching bone screws: $\Phi 6 \times 130\text{mm}$ 4 pcs

Matching instruments: 6mm hex wrench, 6mm screwdriver

S

REF: 6100-0304

Matching bone screws: $\Phi 5 \times 120\text{mm}$ 4 pcs

Matching instruments: 6mm hex wrench, 6mm screwdriver

XS

REF: 6100-0305

Matching bone screws: $\Phi 4 \times 100\text{mm}$ 4 pcs

Matching instruments: 4mm hex wrench, 4mm screwdriver



Dynamic Axial External Fixator(T Type)

XL

REF: 6100-0401

Matching bone screws: $\Phi 6 \times 150\text{mm}$ 2 pcs, HB $\Phi 6.0 \times 150\text{mm}$ 2 pcs

Matching instruments: 6mm hex wrench, 6mm screwdriver

L

REF: 6100-0402

Matching bone screws: $\Phi 6 \times 150\text{mm}$ 2 pcs, HB $\Phi 6.0 \times 150\text{mm}$ 2 pcs

Matching instruments: 6mm hex wrench, 6mm screwdriver

M

REF: 6100-0403

Matching bone screws: $\Phi 6 \times 130\text{mm}$ 2 pcs, HB $\Phi 6.0 \times 150\text{mm}$ 2 pcs

Matching instruments: 6mm hex wrench, 6mm screwdriver

S

REF: 6100-0404

Matching bone screws: $\Phi 5 \times 120\text{mm}$ 2 pcs, HB $\Phi 5.0 \times 130\text{mm}$ 2 pc

Matching instruments: 5mm hex wrench, 5mm screwdriver

XS

REF: 6100-0405

Matching bone screws: $\Phi 4 \times 100\text{mm}$ 4 pcs

Matching instruments: 4mm hex wrench, 4mm screwdriver



Dynamic Axial External Fixator(Ankle Joint Type)

XL

REF: 6100-0501

Matching bone screws: $\Phi 6 \times 150\text{mm}$ 4 pcs

Matching instruments: 6mm hex wrench, 6mm screw driver

L

REF: 6100-0502

Matching bone screws: $\Phi 6 \times 150\text{mm}$ 4 pcs

Matching instruments: 6mm hex wrench, 6mm screw driver

M

REF: 6100-0503

Matching bone screws: $\Phi 6 \times 130\text{mm}$ 4 pcs

Matching instruments: 6mm hex wrench, 6mm screwdriver



Dynamic Axial External Fixator(Hip Joint Type)

XL

REF: 6100-0601

Matching bone screws: $\Phi 6 \times 150\text{mm}$ 2 pcs, HB $\Phi 6.0 \times 150\text{mm}$ 2 pc

Matching instruments: 6mm hex wrench, 6mm screwdriver

L

REF: 6100-0602

Matching bone screws: $\Phi 6 \times 150\text{mm}$ 2 pcs, HB $\Phi 6.0 \times 150\text{mm}$ 2 pc

Matching instruments: 6mm hex wrench, 6mm screwdriver

Pertrochanteric External Fixator

REF: 6100-0701

Matching bone screws: $\Phi 6 \times 150\text{mm}$ 2 pcs

HB $\Phi 6.0 \times 180\text{mm}$ 2 pcs

Matching instruments: 5mm hex wrench, 6mm screwdriver



Elbow Joint External Fixator

REF: 6100-0801

Matching bone screws: $\Phi 5 \times 110\text{mm}$ 4 pcs

Matching instruments: 5mm hex wrench, 3mm hex wrench,
6mm screw driver



Knee Joint External Fixator

REF: 6100-0901

Matching bone screws: $\Phi 6 \times 150\text{mm}$ 2 pcs, $\Phi 6 \times 130\text{mm}$ 2 pcs

Matching instruments: 6mm hex wrench, 3mm hex wrench,
6mm screwdriver



Pelvic External Fixator-Type A

REF: 6100-1001

Matching bone screws: HB $\Phi 6 \times 150\text{mm}$ 4 pcs

Matching instruments: 6mm hex wrench, 6mm screwdriver
Option: 5mm Pin



Pelvic External Fixator-Type B

REF: 6100-1002

Matching bone screws: HB $\Phi 6 \times 150\text{mm}$ 4 pcs

Matching instruments: 6mm hex wrench, 6mm screwdriver
Option: 5mm Pin



Pelvic External Fixator-Type C

REF: 6100-1003

Matching bone screws: HB $\Phi 6 \times 150\text{mm}$ 4 pcs

Matching instruments: 6mm hex wrench, 6mm screwdriver
Option: 5mm Pin



► Bone Lengthening External Fixator-Type A

REF: 6100-1101(250mm) For Children
Matching bone screws: $\Phi 4.0 \times 120\text{mm}$ 6 pcs
Matching instruments: 4mm hex wrench, 4mm screwdriver

REF: 6100-1102(300mm)
Matching bone screws: $\Phi 6 \times 150\text{mm}$ 6 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

REF: 6100-1103(350mm)
Matching bone screws: $\Phi 6 \times 150\text{mm}$ 6 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

REF: 6100-1104(400mm)
Matching bone screws: $\Phi 6 \times 150\text{mm}$ 6 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver



► Bone Lengthening External Fixator-Type B

REF: 6100-1105(250mm)
Matching bone screws: $\Phi 6 \times 150\text{mm}$ 4 pcs (Select 5mm bone screw)
Matching instruments: 6mm hex wrench, 6mm screwdriver

REF: 6100-1106(300mm)
Matching bone screws: $\Phi 6 \times 150\text{mm}$ 4 pcs (Select 5mm bone screw)
Matching instruments: 6mm hex wrench, 6mm screwdriver



► Bone Lengthening External Fixator-Type C

REF: 6100-1107(300mm)
Matching bone screws: $\Phi 6 \times 150\text{mm}$ 6 pcs,
HB $\Phi 6.0 \times 150\text{mm}$ 2 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

► Guide Rail

REF: Z04-01
Size: 200、250、300、400mm



► Intermediate Clamp

REF: Z04-02



► Outer End Clamp

REF: Z04-03



► Adjustment Device

REF: Z04-04
Size: 50、70、110mm



► T-shaped Correction Clamp

REF: Z04-05
Suitable for correcting foot eversion Or varus



► T-clamp

REF: Z04-06



► Elastic Push Block

REF: Z04-07
Apply pressure to the fracture line



► Double Ring Bridge Block

REF: Z04-08
For bridging with ring frame



► Fine Tuning Rotary Clamp

REF: Z04-09
For correction of angular deformity



► Multi-directional Clamp

REF: Z04-10
Suitable for correction of deformed rotation and movement



► Universal T-clamp

REF: Z04-11
For joint end



Transparent Ring Frame Combination (PEEK)

REF: 6100-1402



Transparent Ring Frame Combination(PEEK)

REF: 6100-1401



Tibial Fixator Transparant (PEEK)

REF: 6100-1403



Calcaneal Fixator Transparant (PEEK)

REF: 6100-1404



Ring Connection

REF: HH-01



2/3 Rings

REF: HH-03

Size: Φ125、Φ155、Φ175、
Φ195、Φ215



1/3 Rings


REF: HH-02

Size: Φ125、Φ155、Φ175、
Φ195、Φ215



Foot Rings

REF: HH-04



Rod To Rod Coupling I

REF: HH-05

Size: Φ6/6



Rod To Rod Coupling II

REF: HH-06

Size: Φ6/6



Pin To Ring Coupling

REF: HH-07

Size: Φ1.5/3



Two-hole K-wire Pin

REF: HH-10



90° Extension Plate

REF: HH-11



Single Screw Clamp

REF: HH-12



Diagonal Tie

REF: HH-08



Single Direction Joint With Displacement

REF: HH-13



Ring Bar

REF: HH-09
Size: 50、100、150



Quick Connecting Rod

REF: HH-16



Micro-motion Part For Rings

REF: HH-18



T Wrench

REF: HGJ01



Kirschner Wire Guider

REF: HGJ03



K-wire Pin

REF: HH-14



Single Direction Joint

REF: HH-15



Ring Clamp

REF: HH-17



Kirschner Wire Tensioning Pliers

REF: HGJ02



Open Spanner

REF: HGJ04
Size: #10



Ring External Fixator

REF: 6100-1201
Used for: Tibia, Femur
Configuration:
Half rings 6 pcs, M6×300 threaded rod 3 pcs, M6×20 thread 6 pcs,
M6×20 nut 6 pcs, Pin clamp 16 pcs



2/3 Ring External Fixator

REF: 6100-1202
Used for: Tibia, Femur
Configuration:
Rings with 6 groove 3 pcs, M6×400 threaded
rod 3 pcs, Pin clamp 12 pcs

Lower Limb Lengthening External Fixator

REF: 6100-1203
Used for: Tibia, Femur
Configuration:
Rings 4 pcs, inner threaded connecting rod 4 pcs, M6×180 threaded rod 4
pcs, M6×120 threaded rod 4 pcs, Pin clamp 16 pcs



Ring External Fixator

REF: 6100-1204

Used for: Tibia, Femur

Configuration:

Rings 4 pcs, M6×400 threaded rod 4 pcs, adjusting nut 4 pcs, Pin clamp 16 pcs, Annular tubes 4 pcs



Elbow Joint External Fixator

REF: 6100-1205

Used for: Elbow joint fixator

Configuration:

Rings 4 pcs, Annular diagonal tie 1 pc, single joint 2 pcs, M6×200 threaded rod 8 pcs, Pin clamp 16 pcs, 2-hole connecting plate 1 pc



Ankle External Fixator

REF: 6100-1206

Used for: Equinus

Configuration:

Rings 4 pcs, L connecting 2 pcs, half rings 1 pc, Annular diagonal tie 1 pc, inner threaded post 4 pcs, single joint 2 pcs, M6 threaded rod 8 pcs, Pin clamp 24 pcs

Ring External Fixator

REF: 6100-1207

Used for: Femur, Tibia

Configuration:

Rings with 1 ear, rings with six ears, standard post 6 pcs



Knee Joint External Fixator

REF: 6100-1208

Used for: Knee joint

Configuration:

Rings 4 pcs, foot ring 1 pc, half rings 1 pc, Annular diagonal tie 1 pc, post with 4 edges 2 pcs, single joint 2 pcs, Pin clamp 16 pcs, M6 threaded rod 12 pcs



Ankle Joint External Fixator

REF: 6100-1209

Used for: Equinus

Configuration:

Rings 2 pcs, foot ring 1 pc, L connecting 2 pcs, half rings 1 pc, Annular diagonal tie, Single joint 2 pcs, M6×200 threaded rod 3 pcs, Pin clamp 14 pcs

Parts & Instruments Of Ring External Fixator

Ring With Six Ears

REF: 6100-1210 L
REF: 6100-1211 M
REF: 6100-1212 S
Material: Aluminum



Ring With 3 Ears

REF: 6100-1213 L
REF: 6100-1214 M
REF: 6100-1215 S
Material: Aluminum




Half Ring

REF: HTLL0301 L
REF: HTLL0302 M
REF: HTLL0303 S
Material: Aluminum



Foot Ring

REF: HTLL04 L
Material: Aluminum



Post With 4 Edges

REF: 6100-1216 1 holes
REF: 6100-1217 2 holes
REF: 6100-1218 3 holes
REF: 6100-1219 4 holes
REF: 6100-1220 5 holes
REF: 6100-1221 6 holes
Material: Stainless steel



Post With 6 Edges

REF: 6100-1317
Material: Stainless steel



Rod To Ring Bolt

REF: 6100-1230
Material: Stainless steel



Bush

REF: 6100-1222 ϕ 4
REF: 6100-1223 ϕ 5
REF: 6100-1224 ϕ 6
Material: Aluminum



Pin To Ring Coupling

REF: 6100-1225 ϕ 1.5
REF: 6100-1226 ϕ 2
REF: 6100-1227 ϕ 3
Material: Stainless steel



Pin To Post Coupling

REF: 6100-1228 ϕ 4mm
REF: 6100-1229 ϕ 5mm
Material: Stainless steel



Standard Rod

REF: 6100-1231 L
REF: 6100-1232 M
REF: 6100-1233 S



Packaged Rod


REF: 6100-1234 L
REF: 6100-1235 M
REF: 6100-1236 S



Parts & Instruments Of Ring External Fixator

Rings

REF: 6100-1239 ϕ 100mm
REF: 6100-1240 ϕ 120mm
REF: 6100-1241 ϕ 140mm
REF: 6100-1242 ϕ 160mm
REF: 6100-1243 ϕ 180mm
REF: 6100-1244 ϕ 200mm
REF: 6100-1245 ϕ 220mm
Material: Aluminum



Half Rings

REF: 6100-1247 ϕ 100mm
REF: 6100-1248 ϕ 120mm
REF: 6100-1249 ϕ 140mm
REF: 6100-1250 ϕ 160mm
REF: 6100-1251 ϕ 180mm
REF: 6100-1252 ϕ 200mm
REF: 6100-1253 ϕ 220mm
Material: Aluminum



2/3 Rings

REF: 6100-1255 ϕ 100mm
REF: 6100-1256 ϕ 120mm
REF: 6100-1257 ϕ 140mm
REF: 6100-1258 ϕ 160mm
REF: 6100-1259 ϕ 180mm
REF: 6100-1260 ϕ 200mm
REF: 6100-1261 ϕ 220mm
Material: Aluminum



Foot Ring

REF: 6100-1262 L
REF: 6100-1263 M
REF: 6100-1264 S
Material: Stainless steel



Annular Diagonal Tie

REF: 6100-1237



Threaded Rod M6

REF: 1001-0601/0613(40-400mm)



Adjusting Nut

REF: 6100-1265



Annular Tubes

REF: 6100-1266



Spring

REF: 6100-1267



Single Joint

REF: 6100-1268
Material: Stainless steel



Single Joint

REF: 6100-1269
Material: Stainless steel



Universal Joint

REF: 6100-1270
Material: Stainless steel



Parts & Instruments Of Ring External Fixator

Threaded Column

REF: 6100-1271 1 holes
REF: 6100-1272 2 holes
REF: 6100-1273 3 holes
REF: 6100-1274 4 holes
REF: 6100-1275 5 holes
REF: 6100-1276 6 holes
REF: 6100-1277 7 holes
REF: 6100-1278 8 holes
Material: Stainless steel



Nut Column

REF: 6100-1279 1 holes
REF: 6100-1280 2 holes
REF: 6100-1281 3 holes
REF: 6100-1282 4 holes
REF: 6100-1283 5 holes
REF: 6100-1284 6 holes
REF: 6100-1285 7 holes
REF: 6100-1286 8 holes
Material: Stainless steel



L Connecting Plate

REF: 6100-1287 5 holes
REF: 6100-1288 6 holes
REF: 6100-1289 7 holes
REF: 6100-1290 8 holes
REF: 6100-1291 9 holes
REF: 6100-1292 10 holes
REF: 6100-1293 11 holes
REF: 6100-1294 12 holes
Material: Stainless steel



Connecting Plate

REF: 6100-1295 2 holes
REF: 6100-1296 3 holes
REF: 6100-1297 4 holes
REF: 6100-1298 5 holes
REF: 6100-1299 6 holes
REF: 6100-1300 7 holes
REF: 6100-1301 8 holes
Material: Stainless steel



Cross Connecting Plate

REF: 6100-1302 2 holes
REF: 6100-1303 3 holes
REF: 6100-1304 4 holes
REF: 6100-1305 5 holes
REF: 6100-1306 6 holes
REF: 6100-1307 7 holes
REF: 6100-1308 8 holes
Material: Stainless steel



Inner Threaded Connecting Rod

REF: 6100-1309 80mm
Number: 6100-1310 120mm
Material: Stainless steel



M6 Bolt

REF: HYLL17
Size:
M6×10、12、16、20、25、30
Material: Stainless steel



Pin Clamp

REF: 6100-1311 $\phi 4$
6100-1312 $\phi 5-6$
Material: Stainless steel



Pin Clamp Thick

REF: 6100-1313 $\phi 1.5mm$
REF: 6100-1318 $\phi 2.0mm$
REF: 6100-1314 $\phi 2.5mm$
Material: Stainless steel



Drill Bit Stop Unit

REF: 6100-1315
Material: Stainless steel



Hex Nut M6

REF: 6100-1316
Material: Stainless steel



Wrist Joint External Fixator

REF: HOT001
Used for: wrist
Configuration:
2pcs $\phi 160$ full ring, 3pcs $\phi 160$ half ring, 2pcs hinge, 1pc retractor,
7pcs M6*250 threaded rod



Ankle Joint External Fixator

REF: HOT002
Used for: Ankle joint
Configuration:
4pcs extender, 6pcs M6*150 threaded rod, 3pcs porous plate with
1hole, 3pcs universal joint, 2pcs $\phi 160$ ring, 2pcs $\phi 160U$ ring

Parts & Instruments

Ring

REF: HOTL01
Size:
Φ120、Φ140、Φ150、Φ160、
Φ180、Φ200、Φ220、Φ240
Material: Aluminum




Half Ring

REF: HOTL02
Size:
Φ120、Φ140、Φ150、Φ160、Φ180、
Φ200、Φ220、Φ240
Material: Aluminum



U Ring

REF: HOTL03
Size:
Φ120、Φ140、Φ150、Φ160
Φ180、Φ200
Material: Aluminum



90° Ring

REF: HOTL04
Size:
R90 Long、R90 Short、R120 Long、R120 Short
Material: Aluminum



Porous Plate

REF: HOTL05
Size: 1、2、3、4、5、6、7、8 Holes
Material: Stainless steel



Six Corner Connection

REF: HOTL06
Size: 15mm、30mm
Material: Stainless steel



Connecting Plate

REF: HOTL07
Size: 20、30、40、50
Material: Stainless steel



U Type Ring Extension Rod

REF: HOTL08
Size: 3、4、5、6、7、8 Holes
Material: Aluminum



Ball Head Screw

REF: HOTL09
Size: Φ100、Φ150、Φ200、Φ280
Material: Stainless steel



90° Conversion

REF: HOTL10
Material: Stainless steel



Square Nuts

REF: HOTL11
Material: Stainless steel



Nylon Square Nut

REF: HOTL12
Material: Nylon



Cone Washer

REF: HOTL13
Material: Stainless steel



Cone Washer

REF: HOTL14
Material: Stainless steel



Annular Tubes

REF: HOTL15
Size: 10、20、30、40、50



Parts & Instruments

Connecting Rod

REF: HOTL16
Material: Aluminum



Quick Coupling Rod

REF: HOTL17
Material: Stainless steel



Universal Joint

REF: HOTG03
Material: Stainless steel



Retractor

REF: HOTL18
Size: 100、150、200
Material: Stainless steel



Retractor II

REF: HOTL19
Material: Stainless steel



Threaded Rod

REF: HYLL05



Needle Clamp

REF: HOTJ01
Material: Stainless steel



Square Needle Clamp

REF: HOTJ02
Material: Stainless steel



Pin To Ring Coupling

REF: HOTJ03
Material: Stainless steel



Hinge

REF: HOTG01
Material: Stainless steel



L Type Hinge

REF: HOTG02
Material: Stainless steel



Extender

REF: HOTY01
Size: 70、100、150、200
Material: Aluminum



Tight Wire Clamp

REF: HOTGJ01



Socket Wrench #10

REF: HOTGJ02



Straight Screwdriver #3

REF: HOTGJ03



Screw Driver

REF: HOTGJ04
Size: 6mm



Double Ended Spanner

REF: HOTGJ05
Size: #10



Socket Wrench

REF: HOTGJ06
Size: #10






► Lateral Limb Reconstruction External Fixator

REF: B01
Size: $\phi 4.0 \times 100$, $\phi 4.0 \times 160$

► Large Skin Retractor

REF: Q01
Size: $\phi 2.0 \text{mm} \times 250$
 $\phi 2.5 \text{mm} \times 250$




► Dual-purpose Skin Retractor



REF: Q03

► Small Skin Retractor

REF: Q02
Size: $\phi 1.5 \text{mm} \times 250$
 $\phi 2.0 \text{mm} \times 250$



Hoffmann Type Small Fragment External Fixator Instrument Set 7100-01



Product No.	No.	Product Name	Specification	Qty
7100-0101	1	Rod To Rod Coupling	$\phi 5/5$	12
7100-0102	2	Pin To Rod Coupling	$\phi 5/3-4$	12
7100-0103	3	4-hole Pin Clamp	$\phi 5/3-4$	3
7100-0104	4	Peri-articular Pin Clamp	$\phi 3/4$	1
7100-0105	5	Post	$\phi 5/0^\circ$	2
7100-0106	6	Post	$\phi 5/30^\circ$	4
7100-0107	7	Apex Pin Self-drilling	HA 3.0×80	4
7100-0108			HA 3.0×100	4
7100-0109			HA 4.0×80	4
7100-0110			HA 4.0×100	4
7100-0111			HA 4.0×120	4
7100-0112	8	Drill Sleeves	$\phi 3$	1
7100-0113			$\phi 4$	1
7100-0114	9	Carbon Fiber Connecting Rod	$\phi 5 \times 120$	4
7100-0115			$\phi 5 \times 150$	4
7100-0116			$\phi 5 \times 180$	2
7100-0117			$\phi 5 \times 200$	4
7100-0118	10	T Wrench	#5	1
7100-0119	11	Stabilization/Reduction Wrench	#5	1
7100-0120	12	Elbow Joint Movement Device	$\phi 5$	1
7100-0121	13	Apex Pin Chuck	$\phi 3$	1
7100-0122			$\phi 4$	1
7100-0123	14	T-shaped Pin Wrench	$\phi 3$	1
7100-0124			$\phi 4$	1
7100-0125	15	Hand Drill		1
7100-0126	16	Thumb Wheel	#5/7	1
7100-0127	17	Aluminum Box		1

Parts & Instruments Of Small Fragment External Fixator Instrument Set

Pin To Rod Coupling

REF: 7100-0102
Size: $\Phi 5/3-4$
Material: Aluminum



Rod To Rod Coupling

REF: 7100-0101
Size: $\Phi 5/5$
Material: Aluminum



Pin To Rod Coupling

REF: ZSSP01
Size: $\Phi 5/3-4$
Material: PEEK



Rod To Rod Coupling

REF: ZSSR01
Size: $\Phi 5/5$
Material: PEEK



Pin To Rod Integrated Coupling

REF: ZSSP04
Size: $\Phi 3/4/5$
Material: Aluminum



Peri-articular Pin Clamp

REF: 7100-0104
Size: $\Phi 5/3-4$
Material: Aluminum



4-hole Pin Clamp

REF: 7100-0103
Size: $\Phi 5/3-4$
Material: Aluminum



Pin To Rod Coupling

REF: 7100-0127
Size: $\Phi 3/2$
Material: Aluminum



Rod To Rod Coupling

REF: 7100-0128
Size: $\Phi 3/3$
Material: Aluminum



Parts & Instruments Of Small Fragment External Fixator Instrument Set

Tube To Rod Coupling

REF: 7100-0129
Size: $\Phi 15/5$
Material: Aluminum



Connecting Rod

REF: 7100-0130/7100-0114
Material: Aluminum/Carbon fiber



Straight Post

REF: 7100-0105
Size: $\Phi 5/0^\circ$
Material: Stainless steel



Elbow Joint Movement Device

REF: 7100-0120
Size: $\Phi 5$
Material: Stainless steel



30° post

REF: 7100-0106
Size: $\Phi 5/30^\circ$
Material: Stainless steel



Wrist Joint Movement Device

REF: 7100-0601
Size: $\Phi 5$
Material: Stainless steel



Compression/distraction Tube

REF: 7100-1101
Size: 15mm
Material: Aluminum



Parts & Instruments Of Small Fragment External Fixator Instrument Set

Stabilization/reduction Wrench

REF: 7100-0119
Size: #15



Hand Drill

REF: 7100-0125



Wrench

REF: 7100-0133
Size: #3/#4/#5



T Wrench

REF: 7100-0118
Size: #5



T-shaped Pin Wrench

REF: 7100-0123/7100-0124
Size: Φ3/Φ4



Drill Sleeves

REF: 7100-0112/7100-0113
Size: Φ3、Φ4



Apex Pin Chuck

REF: 7100-0121/7100-0122
Size: Φ3/Φ4



Thumb Wheel

REF: 7100-0126
Size: #5/7




Ulna & Radius External Fixator

REF: 7100-03

Used for: Ulna, Radius

REF	Product Name	Specification	Qty
7100-0301	Rod to rod coupling	Φ5/5	4
7100-0302	4-hole pin clamp	Φ5/3-4	2
7100-0303	Post	Φ5/30°	4
7100-0304	Connecting rod	Φ5×200mm	2
7100-0305	Apex Pin Self-drilling	Φ3×120mm	4



Forearm External Fixator

REF: 7100-04

Used for: Forearm

REF	Product Name	Specification	Qty
7100-0401	Rod to rod coupling	Φ5/5	2
7100-0402	4-hole pin clamp	Φ5/3-4	2
7100-0403	Post	0°	2
7100-0404	Connecting rod	Φ5×200mm	1
7100-0405	Apex Pin Self-drilling	Φ3×120mm	4



Distal Radius External Fixator

REF: 7100-05

Used for: Ulna, Radius

REF	Product Name	Specification	Qty
7100-0501	Rod to rod coupling	Φ5/5	1
7100-0502	Pin to rod coupling	Φ5/3-4	2
7100-0503	Peri-articular pin clamp	Φ5/3-4	1
7100-0504	Connecting rod	Φ5×200mm	1
7100-0505	Apex Pin Self-drilling	Φ3×120mm	4



Wrist Joint External Fixator

REF: 7100-06

Used for: Wrist joint

REF	Product Name	Specification	Qty
7100-0601	Wrist joint movement device	φ5	1
7100-0602	Pin to rod coupling	φ5/3-4	4
7100-0603	Apex Pin Self-drilling	φ3×120mm	4



Elbow Joint External Fixator

REF: 7100-07

Used for: Elbow joint

REF	Product Name	Specification	Qty
7100-0701	Elbow joint movement device	φ5	1
7100-0702	Pin to rod coupling	φ5/3-4	4
7100-0703	Apex Pin Self-drilling	4×120mm	4



Finger & Palm External Fixator

REF: 7100-08

Used for: Finger, Palm

REF	Product Name	Specification	Qty
7100-0801	Pin to rod coupling	φ3/2	2
7100-0802	Connecting rod	φ3×80mm	1
7100-0803	Apex Pin Self-drilling	φ2×70mm	4



Humeral Shaft External Fixator

REF: 7100-09

Used for: Humeral shaft

REF	Product Name	Specification	Qty
7100-0901	Rod to rod coupling	φ5/5	2
7100-0902	4-hole pin clamp	φ5/3-4	2
7100-0903	Post	φ5/30°	2
7100-0904	Connecting rod	φ5×200mm	1
7100-0905	Apex Pin Self-drilling	φ4×120mm	4



Proximal Humerus External Fixator

REF: 7100-10

Used for: Proximal humerus, Condyle neck

REF	Product Name	Specification	Qty
7100-1001	Rod to rod coupling	φ5/5	4
7100-1002	Pin to rod coupling	φ5/3-4	2
7100-1003	4-hole pin clamp	φ5/3-4	1
7100-1004	Connecting rod	φ5×250mm	2
7100-1005	Connecting rod	φ5×200mm	1
7100-1006	Post	φ5/30°	2
7100-1007	Apex Pin Self-drilling	φ4×120mm	4



Humerus External Fixator

REF: 7100-11



Used for: Bone shaft, distal & proximal bone

REF	Product Name	Specification	Qty
7100-1101	Compression/distraction tube	φ15	1
7100-1102	Tube to rod coupling	φ15/5	2
7100-1103	4-hole pin clamp	φ5/3-4	2
7100-1104	Post	φ5/30°	2
7100-1105	Apex Pin Self-drilling	φ4×120mm	4



Hoffmann Type Large Fragment External Fixator Instrument Set 7100-02

31



REF	No.	Product Name	Specification	Qty
7100-0201	1	Rod To Rod Coupling	Φ8/8	12
7100-0202	2	Pin To Rod Coupling	Φ8/5-6	12
7100-0203	3	5-hole Pin Clamp	Φ8/5-6	3
7100-0204	4	Post	Φ8/0°	2
7100-0205	5	Post	Φ8/30°	4
7100-0206			HA 5.0×120	4
7100-0207	6	Apex Pin Self-drilling	HA 5.0×150	8
7100-0208			HA 5.0×180	4
7100-0209	7	Drill Sleeve	Φ5	2
7100-0210			Φ8×200	4
7100-0211	8	Carbon Fiber Connecting Rod	Φ8×250	4
7100-0212			Φ8×300	4
7100-0213	9	Semi-circular Curved Rod	Φ8 180	2
7100-0214			Φ8/ 200	2
7100-0215	10	Stabilization/Reduction Wrench	#7	1
7100-0216	11	Torque Wrench	φ7	1
7100-0217	12	Thumb Wheel	#5/7	1
7100-0218	13	Apex Pin Chuck	φ5	1
7100-0219	14	Hand Drill		1
7100-0220	15	T-shaped Pin Wrench	φ5	1
7100-0221	16	T Wrench	#7	1
7100-0222	17	Aluminum Box		1

Parts & Instruments Of Large Fragment External Fixator Instrument Set

Pin To Rod Coupling

REF: 7100-0201
Size: Φ8/5-6
Material: Aluminum



Rod To Rod Coupling

REF: 7100-0202
Size: Φ8/8
Material: Aluminum



Pin To Rod Coupling

REF: ZSSP01
Size: Φ8/5-6
Material: Aluminum



Rod To Rod Coupling

REF: ZSSR01
Size: Φ8/8
Material: PEEK



Pin To Rod Integrated Coupling

REF: ZSSP04
Size: Φ5/6/8
Material: Aluminum



5-hole Pin Clamp

REF: 7100-0203
Size: Φ8/5-6
Material: Aluminum



Tube To Rod Coupling

REF: 7100-2001
Size: Φ20/8
Material: Aluminum



Semi-circular Curved Rod

REF: 7100-0233/0213/0214
Size: Φ8×160/180/200
Material: Aluminum



Straight Post

REF: 7100-0205
Size: Φ8/0°
Material: Stainless steel



32

HOFFMANN TYPE EXTERNAL FIXATOR SYSTEM

Parts & Instruments Of Large Fragment External Fixator Instrument Set

30° post

REF: 7100-0204
Size: $\Phi 8/30^\circ$
Material: Stainless steel



Connecting Rod

REF: 7100-0234/7100-0210
Material: Aluminum/Carbon fiber



Compression/Distracton Tube

REF: 7100-2002
Size: $\Phi 20$
Material: Aluminum



Compression/Distracton Single Rod

REF: 7100-2028
Size: $\Phi 8$
Material: Aluminum



Compression/Distracton Double Rods

REF: 7100-2029
Size: $\Phi 8$
Material: Aluminum



Ankle Joint Movement Device

REF: 7100-1903
Size: $\Phi 8$
Material: Stainless steel



Knee Joint Movement Device

REF: 7100-1803
Size: $\Phi 8$ L/R
Material: Stainless steel



Parts & Instruments Of Large Fragment External Fixator Instrument Set

Torque Wrench

REF: 7100-0216
Size: #7



Stabilization/Reduction Wrench

REF: 7100-0215
Size: #22



Wrench

REF: 7100-0230
Size: #5/6/7



Hand Drill

REF: 7100-0219



T Wrench

REF: 7100-0220
Size: #7



T-shaped Pin Wrench

REF: 7100-0221
Size: $\Phi 5$, $\Phi 6$



Drill Sleeve

REF: 7100-0209/0231
Size: $\Phi 5/\Phi 6$



Apex Pin Chuck

REF: 7100-0218
Size: $\Phi 5$, $\Phi 6$



Thumb Wheel

REF: 7100-0217
Size: $\Phi 5/\Phi 7$



Proximal Tibia External Fixator

REF: 7100-12

Used for: Proximal tibia

REF	Product Name	Specification	Qty
7100-1201	Rod to rod coupling	φ8/8	4
7100-1202	Pin to rod coupling	φ8/5-6	3
7100-1203	5-hole pin clamp	φ8/5-6	1
7100-1204	Post	φ8/30°	2
7100-1205	Semi-circular Curved Rod	φ8 180	1
7100-1206	Connecting rod	φ8×250mm	2
7100-1207	Apex Pin Self-drilling	φ5×150mm	5




Femoral Shaft External Fixator

REF: 7100-13

Used for: Femoral shaft

REF	Product Name	Specification	Qty
7100-1301	Rod to rod coupling	φ8/8	2
7100-1302	Pin to rod coupling	φ8/5-6	4
7100-1303	Connecting rod	φ8×200mm	1
7100-1304	Connecting rod	φ8×250mm	2
7100-1305	Apex Pin Self-drilling	φ5×150mm	4



Half Ring External Fixator

REF: 7100-14

Used for: Tibia, perone

REF	Product Name	Specification	Qty
7100-1401	Rod to rod coupling	φ8/8	6
7100-1402	5-hole pin clamp	φ8/5-6	2
7100-1403	Post	φ8/30°	4
7100-1404	Connecting rod	φ8×300mm	2
7100-1405	Connecting rod	φ8×250mm	1
7100-1406	Apex Pin Self-drilling	φ5×150mm	4

Proximal Femur External Fixator

REF: 7100-15

Used for: Pertrochanteric

REF	Product Name	Specification	Qty
7100-1501	Rod to rod coupling	φ8/8	4
7100-1502	5-hole pin clamp	φ8/5-6	2
7100-1503	Post	φ8/30°	4
7100-1504	Connecting rod	φ8×300mm	2
7100-1505	Apex Pin Self-drilling	φ5×150mm	2
7100-1506	Apex Pin Self-drilling	φ5×180mm	2




Pelvic External Fixator

REF: 7100-16

Used for: Pelvic

REF	Product Name	Specification	Qty
7100-1601	Rod to rod coupling	φ8/8	6
7100-1602	5-hole pin clamp	φ8/5-6	2
7100-1603	Post	φ8/30°	2
7100-1604	Connecting rod	φ8×300mm	2
7100-1605	Connecting rod	φ8×200mm	2
7100-1606	Apex Pin Self-drilling	φ5×150mm	4



Ankle Joint External Fixator

REF: 7100-17

Used for: Ankle joint

REF	Product Name	Specification	Qty
7100-1701	Rod to rod coupling	φ8/8	3
7100-1702	Pin to rod coupling	φ8/5-6	4
7100-1703	Connecting rod	φ8×250mm	2
7100-1704	Connecting rod	φ8×200mm	1
7100-1705	Apex Pin Self-drilling	φ5×150mm	4



Knee Joint External Fixator

REF: 7100-18

Used for: Knee joint

REF	Product Name	Specification	Qty
7100-1801	Rod to rod coupling	φ8/8	2
7100-1802	Pin to rod coupling	φ8/5-6	4
7100-1803	Knee joint movement device	φ8	1
7100-1804	5-hole pin clamp	φ8/5-6	1
7100-1805	Post	φ8/30°	1
7100-1806	Connecting rod	φ8×150mm	1
7100-1807	Apex Pin Self-drilling	φ5×150mm	6



Ankle Fusion Joint External Fixator

REF: 7100-19

Used for: Ankle joint

REF	Product Name	Specification	Qty
7100-1901	Pin to rod coupling	φ5/3-4	1
7100-1902	Pin to rod coupling	φ8/5-6	4
7100-1903	Ankle joint movement device	φ8	1
7100-1904	Apex Pin Self-drilling	φ5×150mm	4
7100-1905	Apex Pin Self-drilling	φ3×120mm	1



Micro-motion Type External Fixator

REF: 7100-20

REF	Product Name	Specification	Qty
7100-2001	Tube to rod coupling	φ20/8	2
7100-2002	Compression/distraction tube	φ20	1
7100-2003	5-hole pin clamp	φ8/5-6	2
7100-2004	Post	φ8/30°	2
7100-2005	Apex Pin Self-drilling	φ5×150mm	4



Micro-motion Type External Fixator

REF: 7100-21

REF	Product Name	Specification	Qty
7100-2101	Tube-Rod Coupling	φ20/8	2
7100-2102	Compression Distraction Tube	φ20	1
7100-2103	5-Hole Pin Clamp	φ8/5-6	2
7100-2104	Post	φ8/30°	2
7100-2105	Apex Pin Self-drilling	φ6×150mm	4





Synthes Type External Fixator Instrument Set- I 8100-01



REF	No.	Product Name	Specification	Qty
8100-0101	1	Ringent Rod To Rod Coupling	Φ8/8	4
8100-0102	2	Ringent Pin To Rod Coupling	Φ8/4-5	8
8100-0103	3	Ringent Rod To Rod Coupling	Φ11/11	4
8100-0104	4	Ringent Pin To Rod Coupling	Φ11/5-6	12
8100-0105	5	T Clamp	Φ11/5	2
8100-0106	6	Shanz Pin	HA 4.0×90	4
8100-0107			HA 5.0×110	4
8100-0108			HA 5.0×130	6
8100-0109			HA 6.0×130	6
8100-0110			HB 5.0×150	2
8100-0111			HB 6.0×150	2
8100-0112	7	Drill Bits	Φ2.8×200	2
8100-0113			Φ3.8×200	2
8100-0114			Φ4.8×200	2
8100-0115	8	Screw Driver	Φ4	1
8100-0116			Φ5	1
8100-0117	9	Tapered Trocar & Guide	Φ4	1
8100-0118			Φ5	1
8100-0119			Φ6	1
8100-0120	10	Open Wrench	#10	2
8100-0121	11	Pressurizer		1
8100-0423	12	Connecting Rod	Φ8×150	2
8100-0424			Φ8×200	2
8100-0425			Φ8×250	2
8100-0426			Φ8×300	2
8100-0427			Φ11×150	2
8100-0428			Φ11×200	2
8100-0429			Φ11×250	2
8100-0430			Φ11×300	2
8100-0431			Φ11×350	2
8100-0131	13	Aluminum Box		1

Synthes Type External Fixator Instrument Set-II 8100-02



39



REF	No.	Product Name	Specification	Qty
8100-0201	1	Pin To Rod Coupling	Φ11/5-6	12
8100-0202	2	Ringent Rod To Rod Couplin	Φ11/11	9
8100-0203	3	T Clamp	Φ11/5	2
8100-0204	4	Shanz Pin	HA 5.0×110	4
8100-0205			HA 5.0×130	6
8100-0206			HB 5.0×150	2
8100-0207			HA 6.0×130	6
8100-0208			HB 6.0×150	2
8100-0209	5	Drill bits	Φ3.8×200	2
8100-0210			Φ4.8×200	2
8100-0211	6	Screw Driver	Φ5	1
8100-0212			Φ6	1
8100-0213	7	Tapered Trocar & Guide	Φ4	1
8100-0214			Φ5	1
8100-0215			Φ6	1
8100-0216	8	Open Wrench	#10	2
8100-0217	9	Pressurizer		1
8100-0427	10	Connecting Rod	Φ11×150	4
8100-0428			Φ11×200	4
8100-0429			Φ11×250	4
8100-0430			Φ11×300	4
8100-0431			Φ11×350	4
8100-0231	11	Aluminum Box		1

Synthes Type External Fixator Instrument Set-III 8100-03

40



REF	No.	Product Name	Specification	Qty
8100-0301	1	Pin To Rod Coupling	Φ8/3.5-4	8
8100-0302	2	Rod To Rod Coupling	Φ8/8	4
8100-0303	3	Rod To Rod Coupling	Φ11/11	4
8100-0304	4	Pin To Rod Coupling	Φ11/5-6	12
8100-0305	5	T Clamp	Φ11/5	2
8100-0306	6	Shanz Pin	HA 4.0×90	4
8100-0307			HA 5.0×110	4
8100-0308			HA 5.0×130	6
8100-0309			HA 6.0×130	6
8100-0310			HB 5.0×150	2
8100-0311	7	Drill Bits	HB 6.0×150	2
8100-0312			Φ2.8×200	2
8100-0313			Φ3.8×200	2
8100-0314	8	Screw Driver	Φ4.8×200	2
8100-0315			Φ4	1
8100-0316	9	Tapered Trocar & Guide	Φ5	1
8100-0317			Φ4	1
8100-0318			Φ5	1
8100-0319	10	Pressurizer	Φ6	1
8100-0320				1
8100-0321	11	Open Wrench	#10	2
8100-0423	12	Connecting Rod	Φ8×150	2
8100-0424			Φ8×200	2
8100-0425			Φ8×250	2
8100-0426			Φ8×300	2
8100-0427			Φ11×150	2
8100-0428			Φ11×200	2
8100-0429			Φ11×250	2
8100-0430			Φ11×300	2
8100-0431			Φ11×350	2
8100-0331	13	Aluminum Box		1

Parts & Instruments Of Synthes Type External Fixator Instrument Set

Pin To Rod Coupling

REF: 8100-0401
Size: $\Phi 11/5$
Material: Stainless steel



Pin To Rod Coupling

REF: 8100-0402/0403
Size: $\Phi 11/5-6, \Phi 8/3.5-4$
Material: Aluminum



T Clamp

REF: 8100-0411
Size: $\Phi 11/5$
Material: Stainless steel



Rod To Rod Coupling

REF: 8100-0414
Size: $\Phi 11/11$
Material: Stainless steel



Rod To Rod Coupling

REF: 8100-0415/0416
Size: $\Phi 11/11, \Phi 8/8$
Material: Aluminum



T Clamp

REF: 8100-0412/0413
Size: $\Phi 8/4, \Phi 11/5$
Material: Aluminum



Ringent Rod To Rod Coupling

REF: 8100-0404/0405/0406
Size: $\Phi 11/11, \Phi 11/8, \Phi 8/8$
Material: Aluminum



Ringent Pin To Rod Coupling

REF: 8100-0407/0408
Size: $\Phi 8/2.5-3.5, \Phi 8/4-5$
Material: Aluminum



Ringent Pin To Rod Coupling

REF: 8100-0409/0410
Size: $\Phi 11/5-6, \Phi 11/3-4$
Material: Aluminum



Parts & Instruments Of Synthes Type External Fixator Instrument Set

Integrated Pin To Rod Coupling

REF: ZKTP01
Size: $\Phi 8/3-4, \Phi 11/5-6$
Material: Aluminum



Integrated Rod To Rod Coupling

REF: ZKTR01
Size: $\Phi 8/8, \Phi 11/11$
Material: Aluminum



Ringent Rod To Rod Coupling

REF: 8100-0419
Size: $\Phi 11/11$
Material: Aluminum



Pin To Rod Coupling

REF: 8100-0417/0418
Size: $\Phi 8/3-4, \Phi 11/5-6$
Material: Aluminum



Semi-circular Curved Rod

REF: 8100-0420
Size: $\Phi 8\text{mm}$
Material: Aluminum



Trapezoidal Curved Rod

REF: 8100-0422
Size: $\Phi 8\text{mm}$
Material: Aluminum



Connecting Rod

REF: 8100-0423/0431
Size: 8x150/200/250/300
11x150/200/250/300/350



Parts & Instruments Of Synthes Type External Fixator Instrument Set

Curved Rod

REF: 8100-0421
Size: $\Phi 8\text{mm}$
Material: Aluminum



Pressurizer

REF: 8100-0432



Compression/distraction Rod

REF: 8100-0433/0434
Size: $\Phi 8$
Material: Stainless steel



Open Spanner

REF: ZZLGJ01
Size: 10#



Tapered Trocar & Guide

REF: ZZLGJ03
Size: $\Phi 3\text{mm}$, $\Phi 4\text{mm}$, $\Phi 5\text{mm}$, $\Phi 6\text{mm}$
Note: Single or Triple is OK




Hex Wrench

REF: ZZLGJ04
Size: $\Phi 2\text{mm}$, $\Phi 3\text{mm}$, $\Phi 4\text{mm}$, $\Phi 5\text{mm}$, $\Phi 6\text{mm}$



Screw Driver

REF: ZZLGJ05
Size: $\Phi 2\text{mm}$, $\Phi 2.5\text{mm}$, $\Phi 3\text{mm}$, $\Phi 3.5\text{mm}$, $\Phi 4\text{mm}$, $\Phi 5\text{mm}$, $\Phi 6\text{mm}$



Half Ring Type External Fixator

REF: 8100-05

Used for: Tibia & perone open & closed fracture

REF	Product Name	Specification	Qty
8100-0501	Ringent rod to rod coupling	$\Phi 8/8$	6
8100-0502	Ringent pin to rod coupling	$\Phi 8/4-5$	6
8100-0503	Semi-circular curved rod	$\Phi 8$	2
8100-0426	Connecting rod	$\Phi 8 \times 300$	1
8100-0505	Compression/distraction rod	$\Phi 8$	2
8100-0506	Shanz Pin	HA 4.0 $\times$ 120	6



Condylar External Fixator

REF: 8100-06

Used for: Humeral epicondyle, tibial epicondyle, distal fracture

REF	Product Name	Specification	Qty
8100-0601	Ringent rod to rod coupling	$\Phi 8/8$	3
8100-0602	Ringent pin to rod coupling	$\Phi 8/4-5$	6
8100-0603	Semi-circular curved rod	$\Phi 8$	1
8100-0604	Trapezoidal curved rod	$\Phi 8$	1
8100-0425	Connecting rod	$\Phi 8 \times 250$	1
8100-0606	Shanz Pin	HA 4.0 $\times$ 120	6



Tibial Joint Side External Fixator

REF: 8100-07

Used for: Distal tibial joint, proximal fracture

REF	Product Name	Specification	Qty
8100-0701	Ringent rod to rod coupling	$\Phi 8/8$	3
8100-0702	Ringent pin to rod coupling	$\Phi 8/4-5$	7
8100-0703	Semi-circular curved rod	$\Phi 8$	1
8100-0704	Curved rod	$\Phi 8$	1
8100-0426	Connecting rod	$\Phi 8 \times 300$	1
8100-0706	Shanz Pin	HA 4.0 $\times$ 120	7

► Z Frame(T Type) External Fixator

REF: 8100-08

Used for: Femur comminuted fracture

REF	Product Name	Specification	Qty
8100-0801	Ringent rod to rod coupling	φ8/8	2
8100-0802	Ringent pin to rod coupling	φ8/4-5	4
8100-0803	T clamp	φ8/4	1
8100-0804	Connecting rod	φ8×100	1
8100-0423	Connecting rod	φ8×150	2
8100-0806	Shanz Pin	HA 4.0×120	6



► T Frame External Fixator

REF: 8100-09

Used for: Femur, Tibia

REF	Product Name	Specification	Qty
8100-0901	Pin to rod coupling	φ11/5-6	4
8100-0902	T clamp	φ11/5	1
8100-0430	Connecting rod	φ11×300	1
8100-0904	Shanz Pin	φ5×130	4
8100-0905	Shanz Pin	HBφ5×150	2



► Straight Type External Fixator

REF: 8100-10

Used for: Femur, Tibia

REF	Product Name	Specification	Qty
8100-1001	Ringent rod to rod coupling	φ11/11	1
8100-1002	Pin to rod coupling	φ11/5-6	4
8100-0428	Connecting rod	φ11×200	2
8100-1004	Shanz Pin	φ5×130	4



► Z Type External Fixator

REF: 8100-11

Used for: Femur, Tibia

REF	Product Name	Specification	Qty
8100-1101	Rod to rod coupling	φ11/11	2
8100-1102	Pin to rod coupling	φ11/5-6	4
8100-0428	Connecting rod	φ11×200	2
8100-0427	Connecting rod	φ11×150	1
8100-1105	Shanz Pin	φ5×130	4



► T Type External Fixator

REF: 8100-12

Used for: Tibia comminuted fracture

REF	Product Name	Specification	Qty
8100-1201	Pin to rod coupling	φ11/5-6	4
8100-1202	T clamp	φ11/5	1
8100-0430	Connecting rod	φ11×300	1
8100-1204	Shanz Pin	φ5×130	4
8100-1205	Shanz Pin	HBφ5×150	2



► Z Type(Enhancement Mode) External Fixator

REF: 8100-13

Used for: Femur, Tibia

REF	Product Name	Specification	Qty
8100-1301	Rod to rod coupling	φ11/11	2
8100-1302	Pin to rod coupling	φ11/5-6	4
8100-0428	Connecting rod	φ11×200	2
8100-0427	Connecting rod	φ11×150	1
8100-1305	Shanz Pin	φ5×130	4



► Straight Type External Fixator

REF: 8100-14

Used for: Tibia

REF	Product Name	Specification	Qty
8100-1401	Pin to rod coupling	φ11/5-6	4
8100-0429	Connecting rod	φ11×250	1
8100-1403	Shanz Pin	φ5×130	4



► Pelvic Fixator External Fixator

REF: 8100-15

Used for: Pelvic

REF	Product Name	Specification	Qty
8100-1501	Ringent rod to rod coupling	φ11/11	2
8100-1502	Ringent pin to rod coupling	φ11/5-6	4
8100-1503	Ringent rod to rod coupling	φ11/8	2
8100-0428	Connecting rod	φ11×200	2
8100-0430	Connecting rod	φ11×300	1
8100-1506	Compression/distraction rod	φ8	1
8100-1507	Shanz Pin	φ5×150	4



Universal Pin/Rod To Rod Coupling

47

SYNTHESES AO TYPE EXTERNAL FIXATION SYSTEM



› Z Type External Fixator

REF: ZWX001

REF	Product Name	Specification	Qty
ZWXR01	Universal rod to rod coupling	φ11/11	2
ZWXP01	Universal pin to rod coupling	φ11/6	4
ZJYL01	Connecting rod	φ11/200	2
ZJYL01	Connecting rod	φ11×150	1
SS10601345	Threaded rod	φ6×130	4

› Universal Rod To Rod Coupling

REF: ZWXR01
Size: φ11/11、φ8/8



› Universal Pin To Rod Coupling

REF: ZWXP01
Size: φ11/6、φ8/4



Carbon Fiber Composite Ring

48

SYNTHESES AO TYPE EXTERNAL FIXATION SYSTEM




› Three-dimensional Type External Fixator

REF: ZHH001

REF	Product Name	Specification	Qty
ZHHL01	Rings	L	1
ZJYL01	Connecting rod	φ11/250	2
ZHHR01	Rod to ring coupling	φ11/11	2
ZHHP01	Pin to ring coupling	φ2.0	4
ZKKP01	Ringent pin to rod coupling	φ11/5-6	2
ZKKR01	Ringent rod to rod coupling	φ11/11	1
SS10601345	Threaded rod	φ6×130	2
SS5020250	K-wire	φ2.0×250	2

› Tight Wire Clamp

REF: ZHH-G01



› Rings

REF: ZHHL01



› Rod To Ring Coupling

REF: ZHHR01
Size: φ11/11



› Pin To Ring Coupling

REF: ZHHP01
Size: φ2-4



► Straight Type External Fixator

L

REF: 9100-0101

Matching bone screws: $\Phi 6.0 \times 130\text{mm}$ 4 pcs

Matching instruments: 6mm hex wrench, 6mm screw driver

M

REF: 9100-0102

Matching bone screws: $\Phi 6.0 \times 120\text{mm}$ 4 pcs

Matching instruments: 6mm hex wrench, 6mm screw driver
(Select 5mm bone screw)

S

REF: 9100-0103

Matching bone screws: $\Phi 5.0 \times 110\text{mm}$ 4 pcs

Matching instruments: 6mm hex wrench, 5mm screw driver

XS

REF: 9100-0104

Matching bone screws: $\Phi 4.0 \times 90\text{mm}$ 4 pcs

Matching instruments: 4mm hex wrench, 4mm screw driver



► T Type External Fixator

L

REF: 9100-0201

Matching bone screws: $\Phi 6.0 \times 130\text{mm}$ 2 pcs,
HB $\Phi 6.0 \times 150\text{mm}$ 2 pcs

Matching instruments: 6mm hex wrench, 6mm screw driver

M

REF: 9100-0202

Matching bone screws: $\Phi 6.0 \times 120\text{mm}$ 2 pcs,
HB $\Phi 6.0 \times 130\text{mm}$ 2 pcs

Matching instruments: 6mm hex wrench, 6mm screw driver

S

REF: 9100-0203

Matching bone screws: $\Phi 5.0 \times 110\text{mm}$ 2 pcs,
HB $\Phi 5.0 \times 120\text{mm}$ 2 pcs

Matching instruments: 5mm hex wrench, 5mm screw driver

XS

REF: 9100-0204

Matching bone screws: $\Phi 4.0 \times 90\text{mm}$ 4 pcs

Matching instruments: 4mm hex wrench, 4mm screw driver



► Ankle Joint External Fixator

L

REF: 9100-0301

Matching bone screws: $\Phi 6.0 \times 130\text{mm}$ 4 pcs

Matching instruments: 6mm hex wrench, 6mm screw driver

M

REF: 9100-0302

Matching bone screws: $\Phi 5.0 \times 130\text{mm}$ 4 pcs

Matching instruments: 6mm hex wrench, 5mm screw driver

► Straight Type With Middle Pins External Fixator

L

REF: 9100-0401

Matching bone screws: $\Phi 6.0 \times 130\text{mm}$ 6 pcs

Matching instruments: 6mm hex wrench, 6mm screw driver

M

REF: 9100-0402

Matching bone screws: $\Phi 6.0 \times 120\text{mm}$ 6 pcs

Matching instruments: 6mm hex wrench, 6mm screw driver
(Select 5mm bone screw)

S

REF: 9100-0403

Matching bone screws: $\Phi 5.0 \times 110\text{mm}$ 6 pcs

Matching instruments: 6mm hex wrench, 5mm screw driver



► T Type With Middle Pins External Fixator

L

REF: 9100-0501

Matching bone screws: $\Phi 6.0 \times 130\text{mm}$ 4 pcs,
HB $\Phi 6.0 \times 150\text{mm}$ 2 pcs

Matching instruments: 6mm hex wrench, 6mm screw driver

M

REF: 9100-0502

Matching bone screws: $\Phi 6.0 \times 120\text{mm}$ 4 pcs,
HB $\Phi 6.0 \times 130\text{mm}$ 2 pcs

Matching instruments: 6mm hex wrench, 6mm screw driver

S

REF: 9100-0503

Matching bone screws: $\Phi 5.0 \times 110\text{mm}$ 4 pcs,
HB $\Phi 5.0 \times 120\text{mm}$ 2 pcs

Matching instruments: 5mm hex wrench, 5mm screw driver



► 3D Type External Fixator

L

REF: 9100-0601

Matching bone screws: $\Phi 6.0 \times 130\text{mm}$ 4 pcs

Matching instruments: 10-14mm open spanner, 6mm hex wrench

M

REF: 9100-0602

Matching bone screws: $\Phi 5.0 \times 120\text{mm}$ 4 pcs

Matching instruments: 10-14mm open spanner, 6mm hex wrench

S

REF: 9100-0603

Matching bone screws: $\Phi 5.0 \times 110\text{mm}$ 4 pcs

Matching instruments: 10-14mm open spanner, 5mm hex wrench

XS

REF: 9100-0604

Matching bone screws: $\Phi 4.0 \times 90\text{mm}$ 4 pcs

Matching instruments: 10mm open spanner, 4mm hex wrench



Three-section Straight Type External Fixator

REF: 9100-0701

Matching bone screws: $\Phi 6.0 \times 130\text{mm}$ 6 pcs
Matching instruments: 6mm hex wrench,
6mm screw driver



Three-section T-shaped Type External Fixator

REF: 9100-0801

Matching bone screws: $\Phi 6.0 \times 130\text{mm}$ 4 pcs,
HB $\Phi 6.0 \times 150\text{mm}$ 2 pcs
Matching instruments: 6mm hex wrench,
6mm screw driver



Joint Type External Fixator

REF: 9100-0901

Matching bone screws: $\Phi 6.0 \times 130\text{mm}$ 2 pcs,
HB $\Phi 6.0 \times 180\text{mm}$ 2 pcs
Matching instruments: 6mm hex wrench,
6mm screw driver



Oblique Type External Fixator

REF: 9100-1001


Matching bone screws: $\Phi 6.0 \times 150\text{mm}$ 2 pcs,
HB $\Phi 6.0 \times 180\text{mm}$ 2 pcs
Matching instruments: 6mm hex wrench,
6mm screw driver



Dynamic External Fixator

REF: 9100-1101

Matching bone screws: $\Phi 4.0 \times 250\text{mm}$ 2 pcs, HB $\Phi 5.0 \times 250\text{mm}$ 3 pcs
Matching instruments: 10mm open spanner, 4&5 mm screw driver



Stainless Steel Universal Type External Fixator- I

L

REF: 9100-1201

Matching bone screws: $\Phi 6.0 \times 130\text{mm}$ 4 pcs
Matching instruments: 13-16mm open spanner, 4&5 mm hex
wrench, 6mm screw driver

M

REF: 9100-1202

Matching bone screws: $\Phi 5.0 \times 120\text{mm}$ 4 pcs
Matching instruments: 10-14mm open spanner, 4&5 mm hex
wrench, 5mm screw driver

S

REF: 9100-1203

Matching bone screws: $\Phi 5.0 \times 110\text{mm}$ 4 pcs
Matching instruments: 10-14mm open spanner, 4&5 mm hex
wrench, 5 mm screw driver



Stainless Steel Straight Type External Fixator- II

L

REF: 9100-1301

Matching bone screws: $\Phi 6.0 \times 130\text{mm}$ 4 pcs
Matching instruments: 10-16mm open spanner, 4 mm hex
wrench, 6 mm screw driver

M

REF: 9100-1302

Matching bone screws: $\Phi 5.0 \times 120\text{mm}$ 4 pcs
Matching instruments: 10-14mm open spanner, 4 mm hex
wrench, 5 mm screw driver

S

REF: 9100-1303

Matching bone screws: $\Phi 5.0 \times 110\text{mm}$ 4 pcs
Matching instruments: 10-14mm open spanner, 4 mm hex
wrench, 5 mm screw driver



Mini-External Fixator(Finger & Palm)

REF: 9100-1401

Used for: Finger, Palm



Mini- External Fixator(Finger & Palm)

REF: 9100-1402

Used for: Finger, Palm



Wrist Joint External Fixator

REF: 9100-1501

Used for: Wrist joint



Semi-circular Type(Standard) External Fixator

REF: 9100-1601

Used for: Tibia&perone
open&closed fracture

REF	Product Name	Specification	Qty
ZZLR01	Rod to rod coupling	φ8/8	6
ZZLP01	Pin to rod coupling	φ8/4-5	6
ZZLL02	Semi-circular curved rod	φ8mm L	2
ZZLL01	Rod	φ8×300	1
ZZLL04	Compression/distraction rod	φ8	2
SS10401240	Shanz Pin	φ4×120	2
SS5040250	K-wire	φ4×250	2



Pelvic Fixator

REF: 9100-1701

Used for: Pelvic




REF	Product Name	Specification	Qty
ZZLR01	Rod to rod coupling	φ8/8	6
ZZLP01	Pin to rod coupling	φ8/4-5	6
ZZLL03	Semi-circular curved rod	φ8mm	2
ZZLL04	Compression/distraction rod	φ8	1
ZZLL01	Rod	φ8×180	2
SS10501345	Shanz Pin	φ5×130	6

Double Rods Tuberosity External Fixator

REF: 9100-1801

Used for: Pertrochanteric

REF	Product Name	Specification	Qty
ZZLP02	Double rods pin clamp	φ8/4-5	5
ZZLL01	Rod	φ8×300	2
SS10501345	Shanz Pin	φ5×130	2
SS10501850	Shanz Pin	φ5×180	3



Double Rod Single Side External Fixator

REF: 9100-1901

Used for: Femur, Tibia

REF	Product Name	Specification	Qty
ZZLP02	Pin to double rod coupling	φ8/4-5	4
ZZLL01	Rod	φ8×250	2
SS10501345	Shanz Pin	φ5×130	4



Humeral Lateral Condylar Neck External Fixator

REF: 9100-2001

Used for: Humeral lateral condylar neck,
proximal & distal humer

REF	Product Name	Specification	Qty
ZZLP01	Pin to rod coupling	φ8/4-5	4
ZZLR01	Rod to rod coupling	φ8/8	1
ZZLL01	Rod	φ8×200	1
ZZLL01	Rod	φ8×150	1
SS10501345	Shanz Pin	φ5×130	4




Humeral Frame External Fixator

REF: 9100-2101

Used for: Humer

REF	Product Name	Specification	Qty
ZZLR01	Rod to rod coupling	φ8/8	2
ZZLP01	Pin to rod coupling	φ8/4-5	6
ZZLL01	Rod	φ8×250	2
ZZLL01	Rod	φ8×200	1
SS10401240	Shanz Pin	φ4×120	6



► Femoral Condylar External Fixator

REF: 9100-2201

Used for: Femoral supracondylar & intercondylar fracture

REF	Product Name	Specification	Qty
ZZLR01	Rod to rod coupling	φ8/8	3
ZZLP02	Pin to double rod coupling	φ8/4-5	2
ZZLP01	Pin to rod coupling	φ8/4-5	5
ZZLL02	Semi-circular curved rod	φ8mm L	1
ZZLL01	Rod	φ8×300	2
ZZLL03	Curved rod	φ8mm Long	1
SS10501345	Shanz Pin	φ5×130	7



► Femoral Shaft External Fixator

REF: 9100-2301

Used for: Femoral shaft 1/3 to the joint comminuted fracture

REF	Product Name	Specification	Qty
ZZLR01	Rod to rod coupling	φ8/8	3
ZZLP02	Pin to double rod coupling	φ8/4-5	4
ZZLP01	Pin to rod coupling	φ8/4-5	3
ZZLL03	Semi-circular curved rod	φ8mm Long	2
ZZLL01	Rod	φ8×300	2
SS10501345	Shanz Pin	φ5×130	4
SS 5040250	K-wire	φ4×250	1



► Ankle Joint External Fixator

REF: 9100-2401

Used for: Distal tibia, Ankle joint

REF	Product Name	Specification	Qty
ZZLR01	Rod to rod coupling	φ8/8	6
ZZLP01	Pin to rod coupling	φ8-4/5	8
ZZLL04	Compression/distraction rod	φ8	2
ZZLL01	Rod	φ8×350	2
ZZLL01	Rod	φ8×250	2
SS10401240	Shanz Pin	φ4×120	4
SS 5040250	K-wire	φ4×250	2

► Tibial Plateau Frame External Fixator

REF: 9100-2501

Used for: Proximal tibia

REF	Product Name	Specification	Qty
ZZLR01	Rod to rod coupling	φ8/8	3
ZZLP01	Pin to rod coupling	φ8/4-5	8
ZZLL01	Rod	φ8×300	1
ZZLL02	Semi-circular curved rod	φ8mm L	1
ZZLL03	Curved rod	φ8mm Long	1
SS10501345	Shanz Pin	φ5×130	6
SS 5040250	K-wire	φ4×250	1



► Humeral Condylar External Fixator

REF: 9100-2601

Used for: Humeral & tibial lateral epicondyle

REF	Product Name	Specification	Qty
ZZLR01	Rod to rod coupling	φ8/8	3
ZZLP01	Pin to rod coupling	φ8/4-5	5
ZZLL01	Rod	φ8×250	1
ZZLL02	Semi-circular curved rod	φ8mm L	1
ZZLL03	Curved rod	φ8mm Short	1
SS10401240	Shanz Pin	φ4×120	3
SS 5040250	K-wire	φ4×250	1



► Knee Joint External Fixator

REF: 9100-2701

Used for: Spanning joint, distal femur & proximal tibia

REF	Product Name	Specification	Qty
ZZLR01	Rod to rod coupling	φ8/8	8
ZZLP01	Pin to rod coupling	φ8/4-5	14
ZZLL02	Semi-circular curved rod	φ8mm L	2
ZZLL03	Curved rod	φ8mm Long	2
ZZLL06	Spanning joint connecting rod	φ8	2
SS10501345	Shanz Pin	φ5×130	14

Parts And Tools For Combination External Fixator ZZL

Single Rod Clamp

REF: ZZLP01
Size: $\Phi 6/2.5-3.5\text{mm}$ 、 $\Phi 8/3-4\text{mm}$ 、 $\Phi 8/4-5\text{mm}$ 、 $\Phi 8/5-6\text{mm}$



Double Rods Clamp

REF: ZZLR01
Size: $\Phi 6/6\text{mm}$ 、 $\Phi 8/8\text{mm}$



Double Rods Pin Clamp

REF: ZZLP02
Size: $\Phi 8/3-4\text{mm}$ 、 $\Phi 8/4-5\text{mm}$ 、 $\Phi 8/5-6\text{mm}$




Double Rod To Rod Clamp

REF: ZZLR02
Size: $\Phi 6/6\text{mm}$ 、 $\Phi 8/8\text{mm}$



Semi-circular Curved Connecting Rod

REF: ZZLL02
Size: $\Phi 6\text{mm L}$ 、 $\Phi 6\text{mm S}$ 、 $\Phi 8\text{mm L}$ 、 $\Phi 8\text{mm S}$
Material: Stainless steel



Trapezoidal Curved Connecting Rod

REF: ZZLL05
Size: $\Phi 6\text{mm}$ 、 $\Phi 8\text{mm}$
Material: Stainless steel



Curved Connecting Rod

REF: ZZLL03
Size: $\Phi 6\text{mm L}$ 、 $\Phi 6\text{mm S}$ 、 $\Phi 8\text{mm L}$ 、 $\Phi 8\text{mm S}$
Material: Stainless steel



Rod

REF: ZZLL01



Parts And Tools For Combination External Fixator ZZL

Compression/distraction Rod

REF: ZZLL04
Size: $\Phi 8$
Material: Stainless steel



Spanning Joint Connecting Rod

REF: ZZLL06
Size: $\Phi 8$
Material: Stainless steel



Open Spanner

REF: ZZLGJ01
Size: 10#



Open Spanner

REF: ZZLGJ02
Size: 10-13



Tapered Trocar & Guide

REF: ZZLGJ03
Size: $\Phi 3\text{mm}$ 、 $\Phi 4\text{mm}$ 、 $\Phi 5\text{mm}$ 、 $\Phi 6\text{mm}$
Note: Single or Triple is OK




Hex Wrench

REF: ZZLGJ04
Size: $\Phi 2\text{mm}$ 、 $\Phi 3\text{mm}$ 、 $\Phi 4\text{mm}$ 、 $\Phi 5\text{mm}$ 、 $\Phi 6\text{mm}$



Screw Driver

REF: ZZLGJ05
Size: $\Phi 2\text{mm}$ 、 $\Phi 2.5\text{mm}$ 、 $\Phi 3\text{mm}$ 、 $\Phi 3.5\text{mm}$ 、 $\Phi 4\text{mm}$ 、 $\Phi 5\text{mm}$ 、 $\Phi 6\text{mm}$



Hydroxyapatite Coated Shanz Pin



Shanz Pin

REF: 7100-0107/0208
Diameter: 2.0/2.5/3.0/3.5/4.0/5.0/6.0mm
Length: 50-250mm




Apex Pin Self-drilling

Diameter: 3.0/4.0/5.0/6.0mm
Length: 60-250mm



Threaded Guide Wire

Diameter: 1.2/1.5/2.0/2.5/3.0mm
Length: 100/120/150/250/300mm



Wire with Olive

REF: 1001-0401/0402/0403
Specification: 1.5×350/2.0×350/2.5×350mm



K-wire

Diameter: 0.8/1.0/1.2/1.5/1.6/1.8/2.0/2.5/3.0/3.5/4.0/4.5mm
Length: 100/120/150/250/300mm



M6 Threaded Rod

REF: 1001-0601/0613
Length: 40-400mm



Drill Bit

REF: 1001-0701/0734
Diameter: 1.0/1.2/1.5/1.8/2.0/2.5/2.8/3.0/3.2/3.5/4.0/4.2/4.8/5.0mm
Length: 115-200mm



Connecting Rod

Diameter: 3/5/6/8/11mm
Length: 40-400mm
Material: Carbon fiber/Aluminum/Stainless steel

