



Changzhou Meditech Technology Co., Ltd.

Address: 9th. Floor, Inno Laser Building, Science & Education Town, Wujin District,
Changzhou City, Jiangsu Province, China. PC: 213164

Website: www.czmeditech.com

Tel/Fax: +86-519-85855955

HIP & KNEE PROSTHESIS PRODUCT CATALOG

Hip Prosthesis

Bipolar Head Liner
Material: UHMWPE
(Ultrahigh Molecular
Weight Polyethylene)



Description	Code
38/22.225mm	2309.224000438022
40/22.225mm	2309.224000440022
42/28mm	2309.224000442028
44/28mm	2309.224000444028
46/28mm	2309.224000446028
48/28mm	2309.224000448028
50/28mm	2309.224000450028
52/32mm	2309.224000452032
54/32mm	2309.224000454032
56/32mm	2309.224000456032

Bipolar Head Ring
Material: UHMWPE
(Ultrahigh Molecular
Weight Polyethylene)



Description	Code
38/22.225mm	2309.274000038022
40/22.225mm	2309.274000040022
42/28mm	2309.274000042028
44/28mm	2309.274000044028
46/28mm	2309.274000046028
48/28mm	2309.274000048028
50/28mm	2309.274000050028
52/32mm	2309.274000052032
54/32mm	2309.274000054032
56/32mm	2309.274000056032

Bipolar Head Cup
Material: Low Carbon CoCrMo
(Cobalt-Chromium Molybdenum)



Description	Code
38mm	2309.212060038
40mm	2309.212060040
42mm	2309.212060042
44mm	2309.212060044
46mm	2309.212060046
48mm	2309.212060048
50mm	2309.212060050
52mm	2309.212060052
54mm	2309.212060054
56mm	2309.212060056

Hip Prosthesis

Cementless Femoral Stem Short
Material: Ti6AL4V
Surface: HA Coating



Description	Code	Handle length	length neck	OffsetA	neck-shaft angle
1#	241010201	116	30.1	36.4	130
2#	241010202	118	38.4	31.8	130
3#	241010203	121	38.2	31.9	130
4#	241010204	121	39.9	33.4	130
5#	241010205	126	39.7	33.2	130
6#	241010206	128	41.7	35.3	130
7#	241010207	130	41.7	35.3	130
8#	241010208	132	43.7	37.2	130
9#	241010209	134	43.7	37.3	130

Cementless Femoral Stem Long
Material: Ti6AL4V
Surface: HA Coating



Description	Code	Handle length	length neck	OffsetA	neck-shaft angle
0#	2301.241010400	120	28.2	35.3	130
0.5#	2301.2410104005	125	30.1	36.4	130
1#	2301.241010401	130	30.1	36.4	130
2#	2301.241010402	140	31.8	38.4	130
3#	2301.241010403	145	31.9	38.2	130
4#	2301.241010404	150	33.4	39.9	130
5#	2301.241010405	155	33.2	39.7	130
6#	2301.241010406	160	35.2	41.7	130
7#	2301.241010407	165	35.3	41.7	130
8#	2301.241010408	170	37.2	43.7	130

Hip Prosthesis

Cementless Femoral Stem Long
Material: Ti6AL4V
Surface: Titanium Coating



Description	Code	Handle length	length neck	OffsetA	neck-shaft angle
1#	241020201	130	30.1	36.4	130
2#	241020202	140	31.8	38.4	130
3#	241020203	145	31.9	38.2	130
4#	241020204	150	33.4	39.9	130
5#	241020205	155	33.2	39.7	130
6#	241020206	160	35.2	41.7	130
7#	241020207	165	35.3	41.7	130
8#	241020208	170	37.2	43.7	130

Cementless Femoral Stem
Material: Ti6AL4V
Surface: Titanium Coating



Description	Code	Handle length	length neck	OffsetA	neck-shaft angle
2#	2303.241020102	145	32.3	36.7	130
3#	2303.241020103	150	32.3	36.7	130
4#	2303.241020104	155	34.1	38.7	130
5#	2303.241020105	160	34.1	38.7	130
6#	2303.241020106	165	35.9	40.7	130
7#	2303.241020107	170	36	40.7	130
8#	2303.241020108	175	37.8	42.5	130

Hip Prosthesis

Cementless Femoral Revision Stem
Material: Ti6AL4V
Surface: HA Coating



Description	Code	Handle length	length neck	OffsetA	neck-shaft angle
9#	241010509	180	36.1	38.4	130
9.5#	2410105095	185	36.1	38.6	130
10#	241010510	190	37.6	39.95	130
10.5#	2410105105	195	37.6	40	130
11#	241010511	200	39.5	41.7	130

Hip Prosthesis

Metal Femoral Head Material: High Carbon CoCrMo (Cobalt-Chromium-Molybdenum)



Description	Code
22.225+0mm	2307.233060022000
22.225+3mm	2307.233060022003
28-3mm	2307.233060028103
28+0mm	2307.233060028000
28+3mm	2307.233060028003
28+6mm	2307.233060028006
32-4mm	2307.233060032104
32+0mm	2307.233060032000
32+4mm	2307.233060032004
32+7mm	2307.233060032007
36-4mm	2307.233060036104
36+0mm	2307.233060036000
36+4mm	2307.233060036004
36+8mm	2307.233060036008

Metal Femoral Head Material: Low Carbon CoCrMo (Cobalt-Chromium-Molybdenum)



Description	Code
Φ22+0mm	2309.232060022000
Φ22+3mm	2309.232060022003
Φ28-3.5mm	2309.232060028103
Φ28+0mm	2309.232060028000
Φ28+3.5mm	2309.232060028003
Φ28+6mm	2309.232060028006
Φ32-4mm	2309.232060032104
Φ32+0mm	2309.232060032000
Φ32+4mm	2309.232060032004
Φ32+7mm	2309.232060032007

Ceramic Head



Description	Code
28 12/14S	236000028001
28 12/14M	236000028002
28 12/14L	236000028003
32 12/14S	236000032001
32 12/14M	236000032002
32 12/14L	236000032003
32 12/14XL	236000032004
36 12/14S	236000036001
36 12/14M	236000036002
36 12/14L	236000036003
36 12/14XL	236000036004
40 12/14S	236000040001
40 12/14M	236000040002
40 12/14L	236000040003
40 12/14XL	236000040004

Hip Prosthesis

Ceramic Liner



Description	Code
28/35G	226000028035
28/37G	226000028037
32/39G	226000032039
32/41G	226000032041
32/44G	226000032044
32/48G	226000032048
32/52G	226000032052
36/44G	226000036044
36/48G	226000036048
36/52G	226000036052
40/52G	226000040052

XLPE Liner (+4mm 10°)
Material: High Cross Linked Polyethylene



Description	Code
38/22.225	2304.225000338022
40/22.225	2304.225000340022
42/22.225	2304.225000342022
44/22.225	2304.225000344022
46/28	2304.225000346028
48/32	2304.225000348032
50/32	2304.225000350032
52/32	2304.225000352032
54/36	2304.225000354036
56/36	2304.225000356036
58/36	2304.225000358036
60/36	2304.225000360036
62/36	2304.225000362036
64/36	2304.225000364036
66/36	2304.225000366036

PE Liner (+4mm 10°)
Material: UHMWPE
(Ultrahigh Molecular Weight Polyethylene)



Description	Code
38/22	2309.224000338022
40/22	2309.224000340022
42/22	2309.224000342022
44/22	2309.224000344022
46/28	2309.224000346028
48/28	2309.224000348028
50/28	2309.224000350028
52/32	2309.224000352032
54/32	2309.224000354032
56/32	2309.224000356032
58/32	2309.224000358032
60/32	2309.224000360032
62/32	2309.224000362032
64/32	2309.224000364032
66/32	2309.224000366032
68/32	2309.224000368032

Hip Prosthesis

Acetabular Cup-I
Material: Ti6AL4V
Surface: Titanium Coating



Description	Code
38mm	2302.211020138
40mm	2302.211020140
42mm	2302.211020142
44mm	2302.211020144
46mm	2302.211020146
48mm	2302.211020148
50mm	2302.211020150
52mm	2302.211020152
54mm	2302.211020154
56mm	2302.211020156
58mm	2302.211020158
60mm	2302.211020160
62mm	2302.211020162
64mm	2302.211020164
66mm	2302.211020166

Acetabular Cup-II
Material: Ti6AL4V
Surface: Titanium Coating



Description	Code
56mm	2302.211020256
58mm	2302.211020258
60mm	2302.211020260
62mm	2302.211020262
64mm	2302.211020264
66mm	2302.211020266

Acetabular Cup Locking Screw



Description	Code
20×φ6.5mm	2310.261060065020
25×φ6.5mm	2310.261060065025
30×φ6.5mm	2310.261060065030
35×φ6.5mm	2310.261060065035
40×φ6.5mm	2310.261060065040

End Cap



Description	Code
Acetabular Cup End Cap	2310.251060101
Locking Screw End Cap	2310.251060201

Knee Prosthesis

PS Femoral Condyle
Material: CoCrMo
(Cobalt-Chromium-Molybdenum)



Description	Code
2L	2312.10000802
2R	2312.10001002
3L	2312.10000803
3R	2312.10001003
4L	2312.10000804
4R	2312.10001004
5L	2312.10000805
5R	2312.10001005
6L	2312.10000806
6R	2312.10001006
7L	2312.10000807
7R	2312.10001007
8L	2312.10000808
8R	2312.10001008

PS Femoral Insert
Material: UHMWPE
(Ultrahigh Molecular Weight Polyethylene)



Description	Code
2/6	2312.10003502
2/7	2312.10003602
2/8	2312.10003702
2/10	2312.10003802
2/12	2312.10003902
2/14	2312.10004002
2/16	2312.10004102
3/6	2312.10003503
3/7	2312.10003603
3/8	2312.10003703
3/10	2312.10003803
3/12	2312.10003903
3/14	2312.10004003
3/16	2312.10004103

PS Femoral Insert
Material: UHMWPE
(Ultrahigh Molecular Weight Polyethylene)



Description	Code
4/6	2312.10003504
4/7	2312.10003604
4/8	2312.10003704
4/10	2312.10003804
4/12	2312.10003904
4/14	2312.10004004
4/16	2312.10004104
5/6	2312.10003505
5/7	2312.10003605
5/8	2312.10003705
5/10	2312.10003805
5/12	2312.10003905
5/14	2312.10004005
5/16	2312.10004105

Knee Prosthesis

PS Femoral Insert
Material: UHMWPE
(Ultrahigh Molecular
Weight Polyethylene)



Description	Code
6/6	2312.10003506
6/7	2312.10003606
6/8	2312.10003706
6/10	2312.10003806
6/12	2312.10003906
6/14	2312.10004006
6/16	2312.10004106
7/6	2312.10003507
7/7	2312.10003607
7/8	2312.10003707
7/10	2312.10003807
7/12	2312.10003907
7/14	2312.10004007
7/16	2312.10004107

PS Femoral Insert
Material: UHMWPE
(Ultrahigh Molecular
Weight Polyethylene)



Description	Code
8/6	2312.10003508
8/7	2312.10003608
8/8	2312.10003708
8/10	2312.10003808
8/12	2312.10003908
8/14	2312.10004008
8/16	2312.10004108

Knee Prosthesis

CR Femoral Condyle
Material: CoCrMo
(Cobalt-Chromium-Molybdenum)



Description	Code
2L	66.008.1.025
2R	66.008.1.032
3L	66.008.1.026
3R	66.008.1.033
4L	66.008.1.027
4R	66.008.1.034
5L	66.008.1.028
5R	66.008.1.035
6L	66.008.1.029
6R	66.008.1.036

CR Femoral Insert
Material: UHMWPE
(Ultrahigh Molecular Weight Polyethylene)



Description	Code
2/6	66.008.1.000
2/7	66.008.1.001
2/8	66.008.1.002
2/10	66.008.1.003
2/12	66.008.1.004
3/6	66.008.1.005
3/7	66.008.1.006
3/8	66.008.1.007
3/10	66.008.1.008
3/12	66.008.1.009
4/6	66.008.1.010
4/7	66.008.1.011
4/8	66.008.1.012
4/10	66.008.1.013

CR Femoral Insert
Material: UHMWPE
(Ultrahigh Molecular Weight Polyethylene)



Description	Code
4/12	66.008.1.014
5/6	66.008.1.015
5/7	66.008.1.016
5/8	66.008.1.017
5/10	66.008.1.018
5/12	66.008.1.019
6/6	66.008.1.020
6/7	66.008.1.021
6/8	66.008.1.022
6/10	66.008.1.023
6/12	66.008.1.024

Knee Prosthesis

Patella
Material: UHMWPE
(Ultrahigh Molecular
Weight Polyethylene)



Description	Code
2#	2312.10008402
3#	2312.10008403
4#	2312.10008404
5#	2312.10008405

PS/CR Tibial Insert
Material: CoCrMo
(Cobalt-Chromium-Molybdenum)



Description	Code	Size(D/W/H)	Illustration
2L	2312.10008102	44.3 62 37.5	
2R	2312.10008202	44.3 62 37.5	
3L	2312.10008103	46.2 65 37.5	
3R	2312.10008203	46.2 65 37.5	
4L	2312.10008104	48.1 68 39.5	
4R	2312.10008204	48.1 68 39.5	
5L	2312.10008105	50 71 39.5	
5R	2312.10008205	50 71 39.5	
6L	2312.10008106	51.9 74 41.5	
6R	2312.10008206	51.9 74 41.5	
7L	2312.10008107	53.8 77 41.5	
7R	2312.10008207	53.8 77 41.5	
8L	2312.10008108	55.7 80 43.5	
8R	2312.10008208	55.7 80 43.5	

Matching List

Half hip Matching List					
Ceramic Head	Metal Femoral Head(High Carbon)	Metal Femoral Head(Low Carbon)	Bipolar Head Liner	Bipolar Head Ring	Bipolar Head Cup (Low Carbon)
/	22	22	38-22.225mm	38-22.225mm	38mm
/	22	22	40-22.225mm	40-22.225mm	40mm
28	28	28	42-28mm	42-28mm	42mm
28	28	28	44-28mm	44-28mm	44mm
28	28	28	46-28mm	46-28mm	46mm
28	28	28	48-28mm	48-28mm	48mm
28	28	28	50-28mm	50-28mm	50mm
32	32	32	52-32mm	52-32mm	52mm
32	32	32	54-32mm	54-32mm	54mm
32	32	32	56-32mm	56-32mm	56mm

Matching List

Full ceramic Matching List		
Ceramic Head	Ceramic Liner	Acetabular Cup
28	28/35G	44mm
28	28/37G	46mm
32	32/39G	48mm/50mm
32	32/41G	52mm
32	32/44G	54mm/56mm
32	32/48G	58mm/60mm
32	32/52G	62mm/64mm/66mm
36	36/44G	54mm/56mm
36	36/48G	58mm/60mm
36	36/52G	62mm/64mm/66mm
40	40/52G	62mm/64mm/66mm


Half ceramic Matching List		
Ceramic Head	PE Liner	Acetabular Cup
28	46/28	46mm
32	48/32	48mm
32	50/32	50mm
32	52/32	52mm
36	54/36	54mm
36	56/36	56mm
36	58/36	58mm
36	60/36	60mm

Matching List

Metal PE Hip Matching List			
Metal Femoral Head (High Carbon)	Metal Femoral Head (Low Carbon)	PE Liner	Acetabular Cup
22	22	38/22.225	38mm
		40/22.225	40mm
		42/22.225	42mm
		44/22.225	44mm
28	28	46/28	46mm
		48/28	48mm
		50/28	50mm
32	32	48/32	48mm
		50/32	50mm
		52/32	52mm
36	/	54/36	54mm
		56/36	56mm
		58/36	58mm
		60/36	60mm
36	/	62/36	62mm
		64/36	64mm
		66/36	66mm

Instrument Sets

Metal Femoral Head
Instrument Set



Code

DJD2X-2


Acetabular Cup
Instrument Set



Code

DJD2X-1

Cementless Femoral Stem(2303)
Instrument Set




Code

DG91000000


Instrument Sets

Cementless Femoral Stem Short
Instrument Set





Code DG91800000

Ceramic Liner
Instrument Set





Code DG90800000

PS Knee Prosthesis
Instrument Set-I



Code 30.DG93101000

PS Knee Prosthesis
Instrument Set-II



Code 30.DG93101000